




ABHDA Polyclonal Antibody

Catalog No	YP-Ab-05262
Isotype	IgG
Reactivity	Human;Rat
Applications	WB;ELISA
Gene Name	ABHD10
Immunogen	Synthesized peptide derived from human protein . at AA range: 220-300
Specificity	ABHDA Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Abhydrolase domain-containing protein 10;Acyl-protein thioesterase ABHD10; Alpha/beta hydrolase domain-containing protein 10;EC:3.1.1.93;EC:3.1.2.22; ABHDA;ABHD10;mitochondrial (EC 3.4.-.-)
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Melanoma,Ovary,Placenta,
Function	similarity:Belongs to the AB hydrolase superfamily.,
Background	This gene encodes a mitochondrially-localized enzyme that acts in liver cells as a hydrolase. The encoded protein removes glucuronide from mycophenolic acid acyl-glucuronide. There is a pseudogene for this gene on chromosome 6. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night