




ABHDD Polyclonal Antibody

Catalog No	YP-Ab-05263
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ABHD13 C13orf6
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	ABHDD Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABHDD;ABHD13 C13orf6;Abhydrolase domain-containing protein 13 (EC 3.-.-.)
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Brain,Placenta,
Function	sequence caution:Translated as Leu.,similarity:Belongs to the serine esterase family.,
Background	sequence caution:Translated as Leu.,similarity:Belongs to the serine esterase family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night