




ABHD3 Polyclonal Antibody

Catalog No	YP-Ab-05264
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ABHD3 LABH3
Immunogen	Synthesized peptide derived from human protein . at AA range: 240-320
Specificity	ABHD3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ABHD3;ABHD3 LABH3;Abhydrolase domain-containing protein 3 (EC 3.1.1.-) (Lung alpha/beta hydrolase 3)
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Brain,Lung,
Function	similarity:Belongs to the AB hydrolase superfamily. AB hydrolase 4 family.,
Background	This gene encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide range of enzymes. The function of this protein has not been determined. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC12 cells,
primary antibody was diluted at 1:1000, 4° over night