




ACBD4 Polyclonal Antibody

Catalog No	YP-Ab-05272
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ACBD4 HMFT0700
Immunogen	Synthesized peptide derived from human protein . at AA range: 110-190
Specificity	ACBD4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HMFT0700;Acyl-CoA-binding domain-containing protein 4;ACBD4;ACBD4 HMFT0700
Observed Band	29kD
Calculated Molecular Weight	29kD
Cell Pathway	integral component of membrane,
Tissue Specificity	Brain,Colon,Hepatoblastoma,Ovary,Placenta,Thalamus,
Function	similarity:Contains 1 ACB (acyl-CoA-binding) domain.,
Background	This gene encodes a member of the acyl-coenzyme A binding domain containing protein family. All family members contain the conserved acyl-Coenzyme A binding domain, which binds acyl-CoA thiol esters. They are thought to play roles in acyl-CoA dependent lipid metabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night