



AT2L1 Polyclonal Antibody

Catalog No	YP-Ab-05298
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	AGXT2L1
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	AT2L1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AGXT2L1;ETNPPL;AT2L1;Ethanolamine-phosphate phospho-lyase (EC 4.2.3.2) (Alanine--glyoxylate aminotransferase 2-like 1)
Observed Band	50-55 kDa
Calculated Molecular Weight	499 aa, 56 kDa
Cell Pathway	Mitochondrion .
Tissue Specificity	Alzheimer cortex,Brain,Kidney,
Function	cofactor:Pyridoxal phosphate .,similarity:Belongs to the class-III pyridoxal-phosphate-dependent aminotransferase family.,subunit:Homotetramer.,
Background	cofactor:Pyridoxal phosphate .,similarity:Belongs to the class-III pyridoxal-phosphate-dependent aminotransferase family.,subunit:Homotetramer.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
