



CTNL1 Polyclonal Antibody

Catalog No	YP-Ab-05309
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CTNNAL1
Immunogen	Synthesized peptide derived from human protein . at AA range: 150-230
Specificity	CTNL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACRP;Alpha catenin related protein;alpha CATU;Alpha catulin;Alpha-catenin-related protein;CTNL1;CTNNAL1;Alpha-catulin (Alpha-catenin-related protein) (ACRP) (Catenin alpha-like protein 1)
Observed Band	82 kDa
Calculated Molecular Weight	82 kDa
Cell Pathway	Cytoplasm, cytoskeleton . Cell membrane ; Peripheral membrane protein .
Tissue Specificity	Widely expressed. Expressed at lower level in neural tissues and at the highest level in the adrenal gland.
Function	function:May modulate the Rho pathway signaling by providing a scaffold for the Lbc Rho guanine nucleotide exchange factor (ARHGEF1).,induction:Down-regulated in cancer pancreatic cells undergoing differentiation and apoptosis.,similarity:Belongs to the vinculin/alpha-catenin family.,subunit:Interacts with ARHGEF1.,tissue specificity:Widely expressed. Expressed at lower level in neural tissues and at the highest level in the adrenal gland.,
Background	function:May modulate the Rho pathway signaling by providing a scaffold for the Lbc Rho guanine nucleotide exchange factor (ARHGEF1).,induction:Down-regulated in cancer pancreatic cells undergoing differentiation and apoptosis.,



similarity:Belongs to the vinculin/alpha-catenin family.,subunit:Interacts with ARHGEF1.,tissue specificity:Widely expressed. Expressed at lower level in neural tissues and at the highest level in the adrenal gland.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images