



AL2SB Polyclonal Antibody

Catalog No	YP-Ab-05315
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ALS2CR12
Immunogen	Synthesized peptide derived from human protein . at AA range: 250-330
Specificity	AL2SB Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AL2SB;ALS2CR12;Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 12 protein
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Cytoplasm . Cytoplasmic granule . Cell projection, cilium, flagellum . Expressed in the principal piece of the sperm tail, nearest the sperm head. .
Tissue Specificity	Brain,Testis,
Function	regulation of GTPase activity, regulation of hydrolase activity,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
