




ANGL5 Polyclonal Antibody

Catalog No	YP-Ab-05320
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ANGPTL5 UNQ5795/PRO19600
Immunogen	Synthesized peptide derived from human protein . at AA range: 230-310
Specificity	ANGL5 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	angiopoietin like 5;Angiopoietin like protein 5;Angiopoietin related protein 5; ANGPTL5;ANGL5;ANGPTL5 UNQ5795/PRO19600;Angiopoietin-related protein 5 (Angiopoietin-like protein 5)
Observed Band	42kD
Calculated Molecular Weight	42kD
Cell Pathway	Secreted .
Tissue Specificity	Mainly expressed in adult heart.
Function	similarity:Contains 1 fibrinogen C-terminal domain.,tissue specificity:Mainly expressed in adult heart.,
Background	similarity:Contains 1 fibrinogen C-terminal domain.,tissue specificity:Mainly expressed in adult heart.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4°over night