



AQP6 Polyclonal Antibody

Catalog No	YP-Ab-05329
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	AQP6 AQP2L
Immunogen	Synthesized peptide derived from human protein . at AA range: 60-140
Specificity	AQP6 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	240175F9;AQP 6;Aquaporin 6;Aquaporin6;AQP6;AQP6 AQP2L;Aquaporin-6 (AQP-6) (Aquaporin-2-like) (Kidney-specific aquaporin) (hKID)
Observed Band	31kD
Calculated Molecular Weight	31kD
Cell Pathway	Cytoplasmic vesicle membrane ; Multi-pass membrane protein .
Tissue Specificity	Kidney,Testis,
Function	domain:Aquaporins contain two tandem repeats each containing three membrane-spanning domains and a pore-forming loop with the signature motif Asn-Pro-Ala (NPA).,function:Forms a water-specific channel that participates in distinct physiological functions such as glomerular filtration, tubular endocytosis and acid-base metabolism.,similarity:Belongs to the MIP/aquaporin (TC 1.A.8) family.,
Background	The protein encoded by this gene is an aquaporin protein, which functions as a water channel in cells. Aquaporins are a family of small integral membrane proteins related to the major intrinsic protein (MIP or AQP0). This protein is specific for the kidney. This gene and related family members AQP0, AQP2, and AQP5 reside in a cluster on chromosome 12q13. [provided by RefSeq, Jul 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night