



AT7L2 Polyclonal Antibody

Catalog No	YP-Ab-05340
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ATXN7L2
Immunogen	Synthesized peptide derived from human protein . at AA range: 160-240
Specificity	AT7L2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AT7L2;ATXN7L2;Ataxin-7-like protein 2
Observed Band	79kD
Calculated Molecular Weight	79kD
Cell Pathway	
Tissue Specificity	Spleen,
Function	similarity:Contains 1 SCA7 domain.,
Background	similarity:Contains 1 SCA7 domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

