



ABCF3 Polyclonal Antibody

Catalog No	YP-Ab-05354
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	ABCF3
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	ABCF3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATP-binding cassette sub-family F member 3;ABCF3
Observed Band	80 kDa
Calculated Molecular Weight	709 aa, 80 kDa
Cell Pathway	cell-cell adherens junction,membrane,
Tissue Specificity	Cervix,Muscle,Placenta,Skin,Subthalamus nucleus,Vascular sm
Function	caution:Lacks transmembrane domains and is probably not involved in transport.,similarity:Belongs to the ABC transporter family. EF3 subfamily.,similarity:Contains 2 ABC transporter domains.,
Background	caution:Lacks transmembrane domains and is probably not involved in transport.,similarity:Belongs to the ABC transporter family. EF3 subfamily.,similarity:Contains 2 ABC transporter domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images