




BATF2 Polyclonal Antibody

Catalog No	YP-Ab-05369
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	BATF2
Immunogen	Synthesized peptide derived from human protein . at AA range: 180-260
Specificity	BATF2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	B ATF 2;BATF2;Basic leucine zipper transcriptional factor ATF-like 2 (B-ATF-2) (Suppressor of AP-1 regulated by IFN) (SARI)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Nucleus .
Tissue Specificity	Lung,Pericardium,Placenta,Small intestine,
Function	similarity:Belongs to the bZIP family.,similarity:Contains 1 bZIP domain.,
Background	similarity:Belongs to the bZIP family.,similarity:Contains 1 bZIP domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4°over night