



BNIPL Polyclonal Antibody

Catalog No	YP-Ab-05370
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	BNIPL
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	BNIPL Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Bcl-2/adenovirus E1B 19 kDa-interacting protein 2-like protein;BNIP S;BNIPL 1;BNIPL 2;BNIPL1;BNIPL
Observed Band	39kD
Calculated Molecular Weight	39kD
Cell Pathway	nucleus,cytosol,
Tissue Specificity	Isoform 2 is expressed in placenta and lung.
Function	function:May be a bridge molecule between BCL2 and ARHGAP1/CDC42 in promoting cell death.,similarity:Contains 1 CRAL-TRIO domain.,subunit:Homodimer. Interacts with BCL2, ARHGAP1, MIF and GFER.,tissue specificity:Isoform 2 is expressed in placenta and lung.,
Background	The protein encoded by this gene interacts with several other proteins, such as BCL2, ARHGAP1, MIF and GFER. It may function as a bridge molecule between BCL2 and ARHGAP1/CDC42 in promoting cell death. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night