




BRI3 Polyclonal Antibody

Catalog No	YP-Ab-05386
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	BRI3
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	BRI3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BRI3;BRICD2C;Cerebral protein 14;E25;E25C;integral membrane protein 2C;ITM2C;ITM3;PSEC0047;Transmembrane protein BRI3;Brain protein I3 (pRGR2)
Observed Band	13kD
Calculated Molecular Weight	13kD
Cell Pathway	Lysosome membrane ; Multi-pass membrane protein .; [Isoform 1]: Cytoplasm, perinuclear region . Co-localizes with MGAT1 and IFITM3 at the perinuclear region. .; [Isoform 2]: Cytoplasm . Nucleus . Diffuse localization in the cytoplasm and nucleus. .
Tissue Specificity	Brain,Cervix carcinoma,Eye,Gastric adenocarcinoma,Ovary,
Function	similarity:Belongs to the BRI3 family.,
Background	similarity:Belongs to the BRI3 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night