




CAB39 Polyclonal Antibody

Catalog No	YP-Ab-05399
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CAB39 MO25 CGI-66
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CAB39 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein Mo25;MO25;CGI 66;calcium binding protein 39;CAB39;CAB39 MO25 CGI-66;Calcium-binding protein 39 (MO25alpha) (Protein Mo25)
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Duodenum,Hypothalamus,Testis,
Function	function:Together with the STE20-related adaptor-alpha (STRAD alpha) pseudo kinase, forms a regulatory complex capable of stimulating the activity of STK11., similarity:Belongs to the Mo25 family.,
Background	function:Together with the STE20-related adaptor-alpha (STRAD alpha) pseudo kinase, forms a regulatory complex capable of stimulating the activity of STK11., similarity:Belongs to the Mo25 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night