




# CB39L Polyclonal Antibody

|                                    |                                                                                                                                                                                              |
|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-05400                                                                                                                                                                                  |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                          |
| <b>Reactivity</b>                  | Human;Mouse                                                                                                                                                                                  |
| <b>Applications</b>                | WB;ELISA                                                                                                                                                                                     |
| <b>Gene Name</b>                   | CAB39L                                                                                                                                                                                       |
| <b>Immunogen</b>                   | Synthesized peptide derived from part region of human protein                                                                                                                                |
| <b>Specificity</b>                 | CB39L Polyclonal Antibody detects endogenous levels of protein.                                                                                                                              |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                                               |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                       |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                        |
| <b>Dilution</b>                    | WB 1:500-2000 ELISA 1:5000-20000                                                                                                                                                             |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                      |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                         |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                 |
| <b>Synonyms</b>                    | Mo25-like protein;Mo25 like protein;MO25 BETA;Calcium-binding protein 39-like; Antigen MLAA-34;CB39L;CAB39L;Calcium-binding protein 39-like (Antigen MLAA-34) (MO25beta) (Mo25-like protein) |
| <b>Observed Band</b>               | 37kD                                                                                                                                                                                         |
| <b>Calculated Molecular Weight</b> | 37kD                                                                                                                                                                                         |
| <b>Cell Pathway</b>                | cytosol,extracellular exosome,                                                                                                                                                               |
| <b>Tissue Specificity</b>          | Monocytic leukemia,Placenta,                                                                                                                                                                 |
| <b>Function</b>                    | miscellaneous:Found in serum of 50% of patients with acute monocytic leukemia.,similarity:Belongs to the Mo25 family.,                                                                       |
| <b>Background</b>                  | miscellaneous:Found in serum of 50% of patients with acute monocytic leukemia.,similarity:Belongs to the Mo25 family.,                                                                       |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                                                                         |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                            |



## Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night