




CALL6 Polyclonal Antibody

Catalog No	YP-Ab-05409
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	CALML6 CAGLP CALGP
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CALL6 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Calmodulin like protein 6;calmodulin like 6;CALGP;Calglandulin-like protein; Calglandulin like protein;CALL6;CALML6 CAGLP CALGP;Calmodulin-like protein 6 (Calglandulin-like protein)
Observed Band	19kD
Calculated Molecular Weight	19kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Expressed in prostate, thymus, heart, skeleton muscle, bone marrow and ovary.
Function	similarity:Belongs to the calmodulin family. Calglandulin subfamily., similarity:Contains 4 EF-hand domains.,tissue specificity:Expressed in prostate, thymus, heart, skeleton muscle, bone marrow and ovary.,
Background	similarity:Belongs to the calmodulin family. Calglandulin subfamily., similarity:Contains 4 EF-hand domains.,tissue specificity:Expressed in prostate, thymus, heart, skeleton muscle, bone marrow and ovary.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night