



CAN3 Polyclonal Antibody

Catalog No	YP-Ab-05412
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	CAPN3 CANP3 CANPL3 NCL1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CAN3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CAPN3;Calcium-activated neutral proteinase 3;calpain 3;(p94);Calpain L3; Calpain p94;CAN3;CAPN3 CANP3 CANPL3 NCL1;Calpain-3 (EC 3.4.22.54) (Calcium-activated neutral proteinase 3) (CANP 3) (Calpain L3) (Calpain p94) (Muscle-specific calcium-activated neutral protease 3) (New calpain 1) (nCL-1)
Observed Band	9-15 kDa, 55-6 kDa
Calculated Molecular Weight	90kD
Cell Pathway	Cytoplasm. Nucleus, nucleolus .
Tissue Specificity	Isoform I is skeletal muscle specific.
Function	catalytic activity:Broad endopeptidase activity.,disease:Defects in CAPN3 are the cause of limb-girdle muscular dystrophy type 2A (LGMD2A) [MIM:253600]. LGMD2A is an autosomal recessive degenerative myopathy characterized by progressive symmetrical atrophy and weakness of the proximal limb muscles and elevated serum creatine kinase. The symptoms usually begin during the first two decades of life, and the disease gradually worsens, often resulting in loss of walking ability 10 or 20 years after onset.,enzyme regulation:Activated by micromolar concentrations of calcium and inhibited by calpastatin., function:Calcium-regulated non-lysosomal thiol-protease.,online information:Calpain-3 mutations in LGMD2A,similarity:Belongs to the peptidase



C2 family.,similarity:Contains 1 calpain catalytic domain.,similarity:Contains 4 EF-hand domains.,subunit:Interacts with TTN/titin.,tissue specificity:Iso

Background

Calpain, a heterodimer consisting of a large and a small subunit, is a major intracellular protease, although its function has not been well established. This gene encodes a muscle-specific member of the calpain large subunit family that specifically binds to titin. Mutations in this gene are associated with limb-girdle muscular dystrophies type 2A. Alternate promoters and alternative splicing result in multiple transcript variants encoding different isoforms and some variants are ubiquitously expressed. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images