



CEA20 Polyclonal Antibody


Catalog No	YP-Ab-05434
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	CEACAM20 UNQ9366/PRO34155
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CEA20 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CEA20;CEACAM20 UNQ9366/PRO34155;Carcinoembryonic antigen-related cell adhesion molecule 20
Observed Band	64kD
Calculated Molecular Weight	64kD
Cell Pathway	Cell projection, microvillus membrane ; Single-pass type I membrane protein . Apical cell membrane ; Single-pass type I membrane protein . Colocalizes with PTPRH and CEACAM1 at the apical brush border of intestinal cells.
Tissue Specificity	Small intestine,
Function	similarity:Belongs to the immunoglobulin superfamily. CEA family., similarity:Contains 4 Ig-like C2-type (immunoglobulin-like) domains.,
Background	similarity:Belongs to the immunoglobulin superfamily. CEA family., similarity:Contains 4 Ig-like C2-type (immunoglobulin-like) domains.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000