



RABP1 Polyclonal Antibody

Catalog No	YP-Ab-05446
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	CRABP1 RBP5
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	RABP1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRABP;CRABP I;CRABP1;CRABPI;RBP5;RABP1;CRABP1 RBP5;Cellular retinoic acid-binding protein 1 (Cellular retinoic acid-binding protein I) (CRABP-I)
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Brain,
Function	domain:Forms a beta-barrel structure that accommodates hydrophobic ligands in its interior.,function:Cytosolic CRABPs may regulate the access of retinoic acid to the nuclear retinoic acid receptors.,similarity:Belongs to the calycin superfamily. Fatty-acid binding protein (FABP) family.,
Background	This gene encodes a specific binding protein for a vitamin A family member and is thought to play an important role in retinoic acid-mediated differentiation and proliferation processes. It is structurally similar to the cellular retinol-binding proteins, but binds only retinoic acid at specific sites within the nucleus, which may contribute to vitamin A-directed differentiation in epithelial tissue. [provided by RefSeq, Jul 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night