



C1T9B Polyclonal Antibody


Catalog No	YP-Ab-05487
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	C1QTNF9B
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	C1T9B Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1q/TNF-related protein 9B;C1QTNF9B;Complement C1q and tumor necrosis factor-related protein 9B;Complement C1q and tumor necrosis factor-related protein 9-like;CTRP9B;C1T9B;Complement C1q and tumor necrosis factor-related protein 9B (C1q/TNF-related protein 9B) (CTRP9B) (Complement C1q and tumor necrosis factor-related protein 9-like)
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Secreted . Heteromeric complex formation with CTRP9A or ADIPOQ is required for secretion, otherwise, it is retained in the endoplasmic reticulum.
Tissue Specificity	Expressed at low levels. Not expressed in adipose tissues.
Function	similarity:Contains 1 C1q domain.,similarity:Contains 3 collagen-like domains.,
Background	similarity:Contains 1 C1q domain.,similarity:Contains 3 collagen-like domains.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4° over night