



MARH7 Polyclonal Antibody

Catalog No	YP-Ab-05555
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	MARCH7 AXOT RNF177
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	MARH7 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MARH7;MARCH7 AXOT RNF177;E3 ubiquitin-protein ligase MARCH7 (EC 6.3.2.-) (Axotrophin) (Membrane-associated RING finger protein 7) (Membrane-associated RING-CH protein VII) (MARCH-VII) (RING finger protein 177)
Observed Band	77kD
Calculated Molecular Weight	77kD
Cell Pathway	
Tissue Specificity	Liver,Placenta,Skin,Trachea,
Function	domain:The RING-CH-type zinc finger domain is required for E3 ligase activity., function:E3 ubiquitin-protein ligase which may specifically enhance the E2 activity of HIP2. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfer the ubiquitin to targeted substrates.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 RING-CH-type zinc finger.,
Background	MARCH7 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH proteins add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments (Bartee et al., 2004 [PubMed 14722266])



.[supplied by OMIM, Mar 2010],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images