



MARH9 Polyclonal Antibody

Catalog No	YP-Ab-05557
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	MARCH9 RNF179
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	MARH9 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MARH9;MARCH9 RNF179;E3 ubiquitin-protein ligase MARCH9 (EC 6.3.2.-) (Membrane-associated RING finger protein 9) (Membrane-associated RING-CH protein IX) (MARCH-IX) (RING finger protein 179)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein . Lysosome membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitously expressed.
Function	domain:The RING-CH-type zinc finger domain is required for E3 ligase activity., function:E3 ubiquitin-protein ligase that may mediate ubiquitination of MHC-I, CD4 and ICAM1, and promote their subsequent endocytosis and sorting to lysosomes via multivesicular bodies. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfer the ubiquitin to targeted substrates.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 RING-CH-type zinc finger., subunit:Homodimer.,tissue specificity:Ubiquitously expressed.,



Background

MARCH9 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH enzymes add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments. MARCH9 induces internalization of several membrane glycoproteins and directs them to the endosomal compartment (Bartee et al., 2004 [PubMed 14722266]; Hoer et al., 2007 [PubMed 17174307]) [supplied by OMIM, Apr 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night