



ENDD1 Polyclonal Antibody

Catalog No	YP-Ab-05571
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ENDOD1 KIAA0830
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ENDD1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ENDD1;ENDOD1 KIAA0830;Endonuclease domain-containing 1 protein (EC 3.1.30.-)
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Secreted .
Tissue Specificity	Brain,Platelet,
Function	function:May act as a DNase and a RNase .,similarity:Belongs to the DNA/RNA non-specific endonuclease family.,
Background	function:May act as a DNase and a RNase .,similarity:Belongs to the DNA/RNA non-specific endonuclease family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000


