



GCFC1 Polyclonal Antibody

Catalog No	YP-Ab-05598
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	GCFC1 C21orf66 GCFC
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	GCFC1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BM 020;fSAP105;GCFC;GCFC1;GC-rich sequence DNA-binding factor 1; GCFC1 C21orf66 GCFC
Observed Band	93 kDa
Calculated Molecular Weight	917 aa, 105 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous.
Function	function:Possible transcription factor.,similarity:Belongs to the GCF family.,tissue specificity:Ubiquitous.,
Background	PAX3 and PAX7 binding protein 1(PAXBP1) Homo sapiens This gene encodes a protein that may bind to GC-rich DNA sequences, which suggests its involvement in the regulation of transcription. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images