



QSER1 Polyclonal Antibody

Catalog No	YP-Ab-05602
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	QSER1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	QSER1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	glutamine and serine rich 1;QSER1;Glutamine and serine-rich protein 1
Observed Band	190kD
Calculated Molecular Weight	19 kDa
Cell Pathway	Chromosome .
Tissue Specificity	Brain,Epithelium,Testis,Uterus,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
