



GON4L Polyclonal Antibody

Catalog No	YP-Ab-05606
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	GON4L GON4 KIAA1606
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	GON4L Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GON4L;GON4L GON4 KIAA1606;GON-4-like protein (GON-4 homolog)
Observed Band	246kD
Calculated Molecular Weight	246kD
Cell Pathway	Nucleus .
Tissue Specificity	Amygdala,Brain,Cervix carcinoma,Colon,Colon mucosa,Duodenum,Eye,Fetal kidney,Human lung,Liv
Function	function:Enhances transcription activation by YY1.,sequence caution:Translated as Tyr.,similarity:Contains 1 Myb-like domain.,similarity:Contains 2 PAH (paired amphipathic helix) repeats.,subunit:Interacts with YY1.,tissue specificity:Ubiquitous. Detected in small intestine, skeletal muscle, lung, pancreas, brain, stomach, spleen, colon and heart. Detected at very low levels in healthy liver. Highly expressed in most liver carcinomas.,
Background	function:Enhances transcription activation by YY1.,sequence caution:Translated as Tyr.,similarity:Contains 1 Myb-like domain.,similarity:Contains 2 PAH (paired amphipathic helix) repeats.,subunit:Interacts with YY1.,tissue specificity:Ubiquitous. Detected in small intestine, skeletal muscle, lung, pancreas, brain, stomach, spleen, colon and heart. Detected at very low levels



in healthy liver. Highly expressed in most liver carcinomas.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images