



GNA12 Polyclonal Antibody

Catalog No	YP-Ab-05625
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	GNA12
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	GNA12 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	G alpha 12;G alpha-12;G protein subunit alpha 12;gep;GNA 12;GNA12;Guanine nucleotide-binding protein subunit alpha-12 (G alpha-12) (G-protein subunit alpha-12)
Observed Band	41kD
Calculated Molecular Weight	41kD
Cell Pathway	Cell membrane ; Lipid-anchor . Lateral cell membrane ; Lipid-anchor . Cytoplasm . CDH1 enhances cell membrane localization. .
Tissue Specificity	Brain,Pancreas,Placenta,Uterus,
Function	function:Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems., similarity:Belongs to the G-alpha family. G(12) subfamily.,subunit:G proteins are composed of 3 units; alpha, beta and gamma. The alpha chain contains the guanine nucleotide binding site. Interacts with UBXD5.,
Background	function:Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems., similarity:Belongs to the G-alpha family. G(12) subfamily.,subunit:G proteins are composed of 3 units; alpha, beta and gamma. The alpha chain contains the



guanine nucleotide binding site. Interacts with UBXD5.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night