



# GNAO Polyclonal Antibody

|                                    |                                                                                                                                                                                                                                                                                                                                                                                           |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-05629                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Applications</b>                | WB;ELISA                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Gene Name</b>                   | GNAO1                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen</b>                   | Synthesized peptide derived from part region of human protein                                                                                                                                                                                                                                                                                                                             |
| <b>Specificity</b>                 | GNAO Polyclonal Antibody detects endogenous levels of protein.                                                                                                                                                                                                                                                                                                                            |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                                            |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                                     |
| <b>Dilution</b>                    | WB 1:500-2000 ELISA 1:5000-20000                                                                                                                                                                                                                                                                                                                                                          |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Synonyms</b>                    | G ALPHA o;GNAO;GNAO1;Guanine nucleotide-binding protein G(o) subunit alpha                                                                                                                                                                                                                                                                                                                |
| <b>Observed Band</b>               | 38kD                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Calculated Molecular Weight</b> | 38kD                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Cell Pathway</b>                | Cell membrane . Membrane ; Lipid-anchor .                                                                                                                                                                                                                                                                                                                                                 |
| <b>Tissue Specificity</b>          | Brain,Endometrium,Myometrium,                                                                                                                                                                                                                                                                                                                                                             |
| <b>Function</b>                    | function:Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems. The G(o) protein function is not clear.,similarity:Belongs to the G-alpha family. G(i/o/t/z) subfamily.,subunit:G proteins are composed of 3 units; alpha, beta and gamma. The alpha chain contains the guanine nucleotide binding site., |
| <b>Background</b>                  | The protein encoded by this gene represents the alpha subunit of the Go heterotrimeric G-protein signal-transducing complex. Defects in this gene are a cause of early-onset epileptic encephalopathy. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2015]                                                                  |


**matters needing attention**

Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night