




NKX62 Polyclonal Antibody

Catalog No	YP-Ab-05665
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	NKX6-2 GTX NKX6B
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	NKX62 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NKX62;NKX6-2 GTX NKX6B;Homeobox protein Nkx-6.2 (Homeobox protein NK-6 homolog B)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Nucleus .
Tissue Specificity	Highest expression in brain.
Function	similarity:Contains 1 homeobox DNA-binding domain.,tissue specificity:Highest expression in brain.,
Background	similarity:Contains 1 homeobox DNA-binding domain.,tissue specificity:Highest expression in brain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night