



NOBOX Polyclonal Antibody

Catalog No	YP-Ab-05666
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	NOBOX
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	NOBOX Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NOBOX;Homeobox protein NOBOX
Observed Band	76kD
Calculated Molecular Weight	76kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in ovaries, testes and pancreas. Expressed within all stages of the adult female germline, from primordial follicles through to MII oocytes.
Function	
Background	NOBOX oogenesis homeobox(NOBOX) Homo sapiens This homeobox gene encodes a transcription factor that is thought to play a role in oogenesis. In mice, it is essential for folliculogenesis and regulation of oocyte-specific genes. Defects in this gene result in premature ovarian failure type 5.[provided by RefSeq, May 2011],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images