



TF2LY Polyclonal Antibody

Catalog No	YP-Ab-05671
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TGIF2LY TGIFLY
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TF2LY Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TF2LY;TGIF2LY TGIFLY;Homeobox protein TGIF2LY (TGF-beta-induced transcription factor 2-like protein) (TGFB-induced factor 2-like protein;Y-linked) (TGIF-like on the Y)
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Nucleus .
Tissue Specificity	Specifically expressed in adult testis.
Function	function:May have a transcription role in testis. May act as a competitor/regulator of TGIF2LX.,similarity:Belongs to the TALE/TGIF homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,tissue specificity:Specifically expressed in adult testis.,
Background	This gene encodes a member of the TALE/TGIF homeobox family of transcription factors. This gene lies within the male specific region of chromosome Y, in a block of sequence that is thought to be the result of a large X-to-Y transposition. The C-terminus of this protein is divergent from that of its chromosome X homolog (TGIF2LX), suggesting that this protein may act as a regulator of TGIF2LX. [provided by RefSeq, Jul 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night