



UBP54 Polyclonal Antibody

Catalog No	YP-Ab-05679
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	USP54 C10orf29
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	UBP54 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1orf29;Inactive ubiquitin carboxyl-terminal hydrolase 54;Inactive ubiquitin-specific peptidase 54;UBP54;USP54 C10orf29;Inactive ubiquitin carboxyl-terminal hydrolase 54 (Inactive ubiquitin-specific peptidase 54)
Observed Band	185kD
Calculated Molecular Weight	185kD
Cell Pathway	
Tissue Specificity	Weakly expressed in a few tissues.
Function	caution:Although the active site residues are conserved, it lacks the conserved His residue which is normally found 9 residues before the catalytic His., function:Has no peptidase activity.,similarity:Belongs to the peptidase C19 family.,tissue specificity:Weakly expressed in a few tissues.,
Background	caution:Although the active site residues are conserved, it lacks the conserved His residue which is normally found 9 residues before the catalytic His., function:Has no peptidase activity.,similarity:Belongs to the peptidase C19 family.,tissue specificity:Weakly expressed in a few tissues.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000