




LRC42 Polyclonal Antibody

Catalog No	YP-Ab-05690
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	LRRC42
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LRC42 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Leucine-rich repeat-containing protein 42;LRC42;LRRC42
Observed Band	47kD
Calculated Molecular Weight	47kD
Cell Pathway	
Tissue Specificity	Embryo,Placenta,Skin,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the LRRC42 family.,similarity:Contains 2 LRR (leucine-rich) repeats.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the LRRC42 family.,similarity:Contains 2 LRR (leucine-rich) repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4°over night