




TAM41 Polyclonal Antibody

Catalog No	YP-Ab-05740
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TAMM41 C3orf31
Immunogen	Synthesized peptide derived from human protein . at AA range: 220-300
Specificity	TAM41 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C3orf31;DKFZp434E0519;TAM41;TAMM41 C3orf31;Mitochondrial translocator assembly and maintenance protein 41 homolog (TAM41) (MMP37-like protein; mitochondrial)
Observed Band	34kD
Calculated Molecular Weight	34kD
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein ; Matrix side .
Tissue Specificity	Skin,
Function	function:May be involved in the translocation of transit peptide-containing proteins across the mitochondrial inner membrane.,similarity:Belongs to the MMP37 family.,
Background	function:May be involved in the translocation of transit peptide-containing proteins across the mitochondrial inner membrane.,similarity:Belongs to the MMP37 family.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night