



MUC15 Polyclonal Antibody

Catalog No	YP-Ab-05751
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	MUC15 UNQ750/PRO1481
Immunogen	Synthesized peptide derived from human protein . at AA range: 180-260
Specificity	MUC15 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MUC 15;MUC-15;Mucin 15;Mucin-15;PAS3;MUC15;MUC15 UNQ750/PRO1481; Mucin-15 (MUC-15)
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	[Isoform 1]: Cell membrane ; Single-pass type I membrane protein .; [Isoform 2]: Secreted .
Tissue Specificity	Expressed in spleen, thymus, prostate, testis, ovary, small intestine, colon, peripheral blood leukocyte, bone marrow, lymph node and lung.
Function	function:May play a role in the cell adhesion to the extracellular matrix., PTM:Highly glycosylated (N- and O-linked carbohydrates).,sequence caution:Translated as stop.,tissue specificity:Expressed in spleen, thymus, prostate, testis, ovary, small intestine, colon, peripheral blood leukocyte, bone marrow, lymph node and lung.,
Background	function:May play a role in the cell adhesion to the extracellular matrix., PTM:Highly glycosylated (N- and O-linked carbohydrates).,sequence caution:Translated as stop.,tissue specificity:Expressed in spleen, thymus, prostate, testis, ovary, small intestine, colon, peripheral blood leukocyte, bone



marrow, lymph node and lung.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night