




MEGF6 Polyclonal Antibody

Catalog No	YP-Ab-05761
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	MEGF6 EGFL3 KIAA0815
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	MEGF6 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MEGF6;MEGF6 EGFL3 KIAA0815;Multiple epidermal growth factor-like domains protein 6 (Multiple EGF-like domains protein 6) (Epidermal growth factor-like protein 3) (EGF-like protein 3)
Observed Band	169kD
Calculated Molecular Weight	169kD
Cell Pathway	Secreted .
Tissue Specificity	Brain,
Function	similarity:Contains 25 EGF-like domains.,
Background	similarity:Contains 25 EGF-like domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000