




MLF2 Polyclonal Antibody

Catalog No	YP-Ab-05766
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	MLF2
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	MLF2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NTN4;Myeloid leukemia factor 2;Myelodysplasia-myeloid leukemia factor 2;MLF 2;MLF2;Myeloid leukemia factor 2 (Myelodysplasia-myeloid leukemia factor 2)
Observed Band	27kD
Calculated Molecular Weight	27kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Belongs to the MLF family.,tissue specificity:Ubiquitously expressed.,
Background	similarity:Belongs to the MLF family.,tissue specificity:Ubiquitously expressed.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night