



MYO7A Polyclonal Antibody

Catalog No	YP-Ab-05782
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	MYO7A USH1B
Immunogen	Synthesized peptide derived from human protein . at AA range: 830-910
Specificity	MYO7A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DFNA11;DFNB2;myosin VIIA;MYOVIIA;MYU7A;MYO7A;MYO7A USH1B; Unconventional myosin-VIIa
Observed Band	240-250 kDa
Calculated Molecular Weight	254 kDa
Cell Pathway	Cytoplasm . Cytoplasm, cell cortex . Cytoplasm, cytoskeleton . Cell junction, synapse . In the photoreceptor cells, mainly localized in the inner and base of outer segments as well as in the synaptic ending region (PubMed:8842737). In retinal pigment epithelial cells colocalizes with a subset of melanosomes, displays predominant localization to stress fiber-like structures and some localization to cytoplasmic puncta (PubMed:19643958, PubMed:27331610). Detected at the tip of cochlear hair cell stereocilia (PubMed:21709241). The complex formed by MYO7A, USH1C and USH1G colocalizes with F-actin (PubMed:21709241). .
Tissue Specificity	Expressed in the pigment epithelium and the photoreceptor cells of the retina. Also found in kidney, liver, testis, cochlea, lymphocytes. Not expressed in brain.
Function	alternative products:Additional isoforms seem to exist,developmental stage:Detected in optic cup in 5.5 weeks-old embryos. Expressed in retinal pigment epithelium, cochlear and vestibular neuroepithelia, and olfactory



epithelium at 8 weeks. At 19 weeks, present in both pigment epithelium and photoreceptor cells. At 24-28 weeks, expression in pigment epithelium and photoreceptor cells increases. Present in pigment epithelium and photoreceptor cells in adult.,disease:Defects in MYO7A are the cause of non-syndromic sensorineural deafness autosomal dominant type 11 (DFNA11) [MIM:601317]., disease:Defects in MYO7A are the cause of non-syndromic sensorineural deafness autosomal recessive type 2 (DFNB2) [MIM:600060]; also called neurosensory non-syndromic recessive deafness 2 (NSRD2). DFNB2 is a form of sensorineural hearing loss. Sensorineural deafness results from damage to the neural receptor

Background

This gene is a member of the myosin gene family. Myosins are mechanochemical proteins characterized by the presence of a motor domain, an actin-binding domain, a neck domain that interacts with other proteins, and a tail domain that serves as an anchor. This gene encodes an unconventional myosin with a very short tail. Defects in this gene are associated with the mouse shaker-1 phenotype and the human Usher syndrome 1B which are characterized by deafness, reduced vestibular function, and (in human) retinal degeneration. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images