



MYO15 Polyclonal Antibody

Catalog No	YP-Ab-05784
Isotype	IgG
Reactivity	Human;Mouse
Applications	IHC;IF
Gene Name	MYO15A MYO15
Immunogen	Synthesized peptide derived from human protein . at AA range: 2990-3070
Specificity	MYO15 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC-p 1:50-300. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MYO15;MYO15A MYO15;Unconventional myosin-XV (Unconventional myosin-15)
Observed Band	388kD
Calculated Molecular Weight	388kD
Cell Pathway	Cell projection, stereocilium . Cytoplasm, cytoskeleton . Localizes to stereocilium tips in cochlear and vestibular hair cells. .
Tissue Specificity	Highly expressed in pituitary. Also expressed at lower levels in adult brain, kidney, liver, lung, pancreas, placenta and skeletal muscle. Not expressed in brain. In the pituitary, highly expressed in anterior gland cells.
Function	disease:Defects in MYO15A are the cause of non-syndromic sensorineural deafness autosomal recessive type 3 (DFNB3) [MIM:600316]. DFNB3 is a form of sensorineural hearing loss. Sensorineural deafness results from damage to the neural receptors of the inner ear, the nerve pathways to the brain, or the area of the brain that receives sound information.,function:Myosins are actin-based motor molecules with ATPase activity. Unconventional myosins serve in intracellular movements. Their highly divergent tails are presumed to bind to membranous compartments, which would be moved relative to actin filaments. Required for the arrangement of stereocilia in mature hair bundles.,



similarity:Contains 1 FERM domain.,similarity:Contains 1 myosin head-like domain.,similarity:Contains 1 SH3 domain.,similarity:Contains 2 MyTH4 domains.,similarity:Contains 3 IQ domains.,subcellular location:Localizes to ste

Background

This gene encodes an unconventional myosin. This protein differs from other myosins in that it has a long N-terminal extension preceding the conserved motor domain. Studies in mice suggest that this protein is necessary for actin organization in the hair cells of the cochlea. Mutations in this gene have been associated with profound, congenital, neurosensory, nonsyndromal deafness. This gene is located within the Smith-Magenis syndrome region on chromosome 17. Read-through transcripts containing an upstream gene and this gene have been identified, but they are not thought to encode a fusion protein. Several alternatively spliced transcript variants have been described, but their full length sequences have not been determined. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded human spleen tissue. 1,primary Antibody was diluted at 1:200(4° overnight). 2, Sodium citrate pH 6.0 was used for antigen retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200