



NRN1 Polyclonal Antibody

Catalog No	YP-Ab-05827
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	NRN1 NRN
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	NRN1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NRN1;NRN1 NRN;Neuritin
Observed Band	15kD
Calculated Molecular Weight	15kD
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor . Cell junction, synapse .
Tissue Specificity	Adrenal cortex,Brain,Uterus,
Function	function:Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.,similarity:Belongs to the neuritin family.,
Background	This gene encodes a member of the neuritin family, and is expressed in postmitotic-differentiating neurons of the developmental nervous system and neuronal structures associated with plasticity in the adult. The expression of this gene can be induced by neural activity and neurotrophins. The encoded protein contains a consensus cleavage signal found in glycosylphosphatidylinositol (GPI)-anchored proteins. The encoded protein promotes neurite outgrowth and arborization, suggesting its role in promoting neuritogenesis. Overexpression of the encoded protein may be associated with astrocytoma progression. Alternative splicing results in multiple transcript variants. [provided by RefSeq,




Jul 2013],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from Jarkat cells,
primary antibody was diluted at 1:1000, 4°over night