




NPIP Polyclonal Antibody

Catalog No	YP-Ab-05868
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	NPIP
Immunogen	Synthesized peptide derived from human protein . at AA range: 110-190
Specificity	NPIP Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NPIP;Nuclear pore complex-interacting protein (NPIP)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Nucleus, nuclear pore complex . Nucleus membrane . Colocalizes with nuclear pore complex protein NUP62.
Tissue Specificity	Widely expressed.
Function	similarity:Belongs to the NPIP family.,similarity:Contains 1 GPS domain., subcellular location:Colocalizes with nuclear pore complex protein NUP62., subunit:May associate with the nuclear pore complex.,tissue specificity:Widely expressed.,
Background	similarity:Belongs to the NPIP family.,similarity:Contains 1 GPS domain., subcellular location:Colocalizes with nuclear pore complex protein NUP62., subunit:May associate with the nuclear pore complex.,tissue specificity:Widely expressed.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night.