



P210L Polyclonal Antibody

Catalog No	YP-Ab-05870
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	NUP210L
Immunogen	Synthesized peptide derived from human protein . at AA range: 980-1060
Specificity	P210L Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	P210L;NUP210L;Nuclear pore membrane glycoprotein 210-like (Nucleoporin 210 kDa-like) (Nucleoporin Nup210-like)
Observed Band	207kD
Calculated Molecular Weight	207kD
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Testis,
Function	similarity:Belongs to the NUP210 family.,
Background	similarity:Belongs to the NUP210 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
