



NRBP2 Polyclonal Antibody

Catalog No	YP-Ab-05874
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	NRBP2 PP9320 TRG16
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	NRBP2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NRBP2;pp9320;TRG 16;TRG16;NRBP2 PP9320 TRG16;Nuclear receptor-binding protein 2 (Transformation-related gene 16 protein) (TRG-16)
Observed Band	33 kDa, 58 kDa
Calculated Molecular Weight	501 aa, 58 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Testis,
Function	domain:The protein kinase domain is predicted to be catalytically inactive., sequence caution:Intron retention.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,
Background	domain:The protein kinase domain is predicted to be catalytically inactive., sequence caution:Intron retention.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images