



NOC4L Polyclonal Antibody

Catalog No	YP-Ab-05877
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	NOC4L
Immunogen	Synthesized peptide derived from human protein . at AA range: 370-450
Specificity	NOC4L Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NET49;NOC4;NOC4-like protein;Nucleolar complex-associated protein 4-like protein;NOC4L;Nucleolar complex protein 4 homolog (NOC4 protein homolog) (NOC4-like protein) (Nucleolar complex-associated protein 4-like protein)
Observed Band	59 kDa
Calculated Molecular Weight	516 aa, 59 kDa
Cell Pathway	Nucleus membrane ; Multi-pass membrane protein . Nucleus, nucleolus .
Tissue Specificity	Brain, Eye, Ovary,
Function	similarity: Belongs to the CBF/MAK21 family.,
Background	similarity: Belongs to the CBF/MAK21 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
