



NOG1 Polyclonal Antibody

Catalog No	YP-Ab-05878
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	GTPBP4 CRFG NOG1
Immunogen	Synthesized peptide derived from human protein . at AA range: 580-660
Specificity	NOG1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NGB;GTP-binding protein NGB;GTP-binding protein 4;GTP binding protein NGB;GTP binding protein 4;NOG1;GTPBP4 CRFG NOG1;Nucleolar GTP-binding protein 1 (Chronic renal failure gene protein) (GTP-binding protein NGB)
Observed Band	74 kDa
Calculated Molecular Weight	634 aa, 74 kDa
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Cervix carcinoma,Epithelium,Kidney,Liver,Uterus,
Function	function:Involved in the biogenesis of the 60S ribosomal subunit., similarity:Belongs to the GTP1/OBG family. NOG subfamily.,
Background	GTP-binding proteins are GTPases and function as molecular switches that can flip between two states: active, when GTP is bound, and inactive, when GDP is bound. ‘Active’ in this context usually means that the molecule acts as a signal to trigger other events in the cell. When an extracellular ligand binds to a G-protein-linked receptor, the receptor changes its conformation and switches on the trimeric G proteins that associate with it by causing them to eject their GDP and replace it with GTP. The switch is turned off when the G protein hydrolyzes its own bound GTP, converting it back to GDP. But before

that occurs, the active protein has an opportunity to diffuse away from the receptor and deliver its message for a prolonged period to its downstream target. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images