



NOL11 Polyclonal Antibody

Catalog No	YP-Ab-05880
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	NOL11 L14
Immunogen	Synthesized peptide derived from human protein . at AA range: 510-590
Specificity	NOL11 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	L14;NOL11;nucleolar protein 11;NOL11 L14
Observed Band	75-81 kDa
Calculated Molecular Weight	719 aa, 81 kDa
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Lung,Placenta,Uterus,
Function	
Background	NOL11 (Nucleolar Protein 11) is a Protein Coding gene. Among its related pathways are rRNA processing in the nucleus and cytosol and Gene Expression. Ribosome biogenesis factor. May be required for both optimal rDNA transcription and small subunit (SSU) pre-rRNA processing at sites A, A0, 1 and 2b.
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images