



NOL9 Polyclonal Antibody

Catalog No	YP-Ab-05885
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	NOL9
Immunogen	Synthesized peptide derived from human protein . at AA range: 640-720
Specificity	NOL9 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NOL9;nucleolar protein 9;Polynucleotide 5'-hydroxyl-kinase NOL9 (EC 2.7.1.-) (Nucleolar protein 9)
Observed Band	79 kDa
Calculated Molecular Weight	702 aa, 79 kDa
Cell Pathway	Nucleus . Nucleus, nucleolus . Colocalizes with pre-60S rRNP particles.
Tissue Specificity	Amygdala,Hepatoma,Kidney,Lymph,Ovary,Skin,
Function	similarity:Belongs to the NOL9 family.,
Background	similarity:Belongs to the NOL9 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images
