



PHX2A Polyclonal Antibody

Catalog No	YP-Ab-05897
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	PHOX2A ARIX PMX2A
Immunogen	Synthesized peptide derived from human protein . at AA range: 210-290
Specificity	PHX2A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARIX;ARIX1 homeodomain protein;CFEOM2;FEOM2;paired like homeobox 2a; PHOX2A;PMX2A;PHX2A;PHOX2A ARIX PMX2A;Paired mesoderm homeobox protein 2A (ARIX1 homeodomain protein) (Aristaless homeobox protein homolog) (Paired-like homeobox 2A)
Observed Band	31kD
Calculated Molecular Weight	31kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,
Function	disease:Defects in PHOX2A are the cause of congenital fibrosis of extraocular muscles type 2 (CFEOM2) [MIM:602078]. CFEOM encompasses several different inherited strabismus syndromes characterized by congenital restrictive ophthalmoplegia affecting extraocular muscles innervated by the oculomotor and/or trochlear nerves. CFEOM is characterized clinically by anchoring of the eyes in downward gaze, ptosis, and backward tilt of the head. CFEOM2 may result from the aberrant development of the oculomotor (nIII), trochlear (nIV) and abducens (nVI) cranial nerve nuclei.,function:May be involved in regulating the specificity of expression of the catecholamine biosynthetic genes. Acts as a transcription activator/factor. Could maintain the noradrenergic phenotype.,



similarity:Belongs to the paired homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,

Background

paired like homeobox 2a(PHOX2A) Homo sapiens The protein encoded by this gene contains a paired-like homeodomain most similar to that of the Drosophila aristaless gene product. The encoded protein plays a central role in development of the autonomic nervous system. It regulates the expression of tyrosine hydroxylase and dopamine beta-hydroxylase, two catecholaminergic biosynthetic enzymes essential for the differentiation and maintenance of the noradrenergic neurotransmitter phenotype. The encoded protein has also been shown to regulate transcription of the alpha3 nicotinic acetylcholine receptor gene. Mutations in this gene have been associated with autosomal recessive congenital fibrosis of the extraocular muscles. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night