



TYSD1 Polyclonal Antibody

Catalog No	YP-Ab-05904
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TYSND1
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	TYSD1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TYSD1;TYSND1;Peroxisomal leader peptide-processing protease (EC 3.4.21.-) (Trypsin domain-containing protein 1) [Cleaved into: Peroxisomal leader peptide-processing protease;15 kDa form;Peroxisomal leader peptid
Observed Band	62kD
Calculated Molecular Weight	62kD
Cell Pathway	Peroxisome .
Tissue Specificity	Brain,Lung,
Function	function:Peroxisomal protease that mediates both the removal of the leader peptide from proteins containing a PTS2 target sequence and the specific processing of PTS1 proteins. May participate to a regulatory mechanism that control the peroxisomal beta-oxidation of fatty acids.,PTM:Processed into the 49 kDa and the 10 kDa forms upon import into the peroxisomes. The 49 kDa form probably corresponds to the mature enzyme, while the 10 kDa form may represent an inhibitory N-terminal fragment which prevents protease activity until it reaches the peroxisomes.,similarity:Belongs to the peptidase S1B family.,

**Background**

This gene encodes a protease that removes the N-terminal peroxisomal targeting signal (PTS2) from proteins produced in the cytosol, thereby facilitating their import into the peroxisome. The encoded protein is also capable of removing the C-terminal peroxisomal targeting signal (PTS1) from proteins in the peroxisomal matrix. The full-length protein undergoes self-cleavage to produce shorter, potentially inactive, peptides. Alternative splicing results in multiple transcript variants for this gene. [provided by RefSeq, Jan 2013],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images